

## SET-UP SAMPLE

### IMZ S 501

#### HIGH ALLOY STEEL

Sample size: approx. fi 47 mm x 23 mm

Chemical composition, % m/m:

|    |        |
|----|--------|
| C  | 0.15   |
| Mn | 0.32   |
| Si | 0.40   |
| P  | 0.02   |
| S  | <0.005 |
| Cr | 13.2   |
| Ni | 0.11   |

Approved by  
Director of the Institute

Prof. Dr. Hab. Eng. Adam Zieliński

Issue date: 02 July 2025

Revision history:

02 July 2025 (editorial changes)

December 2021 (Original document date)